

# epinephrine<sup>®</sup> snap-EMS

EPINEPHRINE CONVENIENCE KIT



## Each kit contains:

- One Vial of Epinephrine 1 mg/mL
- Two Epinephrine Syringes™
- Two 23-Gauge, 1-Inch Safety Needles
- Four Alcohol Prep Pads
- One Epinephrine Package Insert
- One Instructions For Use Card

Epinephrine Snap<sup>®</sup>-EMS is a complete, ready-to-use epinephrine convenience kit engineered for a fast, accurate response to anaphylactic emergencies.

Each kit contains one vial of Adrenalin 1 mg/mL and all essential supplies for safe, precise administration during an emergency.

Epinephrine Snap<sup>®</sup>-EMS is designed to treat anaphylaxis in patients weighing 15 kg (33 lbs) or greater.

Protected by U.S. Patent No. 11,833,336 and 12,589,210, featuring our innovative syringe and dosing system engineered to deliver predetermined, accurate epinephrine doses—improving safety, speed, compliance, and cost efficiency when every second matters.

FDA Registered. Made in the USA.



## The Epinephrine Syringe™

Our proprietary syringe features clear, easy-to-read indicators at 0.15 mL, 0.3 mL, and 0.5 mL, ensuring fast and accurate dose selection. When used with Epinephrine 1 mg/mL [included in every Epinephrine Snap<sup>®</sup>-EMS kit] the syringe enables precise administration of 0.15 mg, 0.3 mg, and 0.5 mg epinephrine doses during anaphylactic emergencies.

Learn more at [snapmedicalindustries.com](http://snapmedicalindustries.com).

The FIRST Epinephrine Convenience Kit  
Available for Nationwide Delivery

Low Cost. High Value. Always Prepared.