

the epinephrine syringe™



Description:

- For Professional use by health care providers.
- Luer-lock syringe intended for use with epinephrine 1 mg/mL only for intramuscular or subcutaneous injection. Not for intravenous use.
- Features dose markings at 0.15 mL, 0.3 mL, and 0.5 mL.
- Indicated for use in patients weighing ≥ 15 kg (33 lbs).
- Single-use only. Discard after use in an appropriate container.
- Manufactured without natural rubber latex. Non-toxic and pyrogen-free.
- Supplied without needle.
- Available in boxes of 100 individually packaged syringes.

Precision you can trust—when it matters most.

Always follow your anaphylaxis treatment order/protocol when deciding if The Epinephrine Syringe™ meets your needs in treating anaphylaxis.

Learn more at snapmedicalindustries.com.

Protected by US Patent number 11833336

FDA Registered

